

SPRI Bead Purification Protocol

SPRI beads are Agencourt AMPure XP beads (Cat. No. A63881).

- 1** SPRI beads are stored at 4 C. Remove from refrigerator and warm to room temperature.
- 2** Transfer the reaction to a 1.7 mL tube, because that is what fits on the magnet.
- 3** Calculate volume of beads needed: $1.8 \times \text{rxn volume}$. I.e. for 50 uL rxn, use 90 uL beads.
- 4** Measure that volume in a pipette tip with water so you know the correct level.
- 5** Mix beads thoroughly by vortexing. Then pipette the needed volume into the reaction.
- 6** Mix the beads and the reaction by flicking and gently vortexing. Then quickly spin any solution down to the bottom of the tube (only 1 or 2 seconds -- avoid pelleting the beads).
- 7** Let stand at room temperature for 5 minutes.
During this time make $1.1 \text{ mL} \times \# \text{ rxns}$ of 80% EtOH.
- 8** Place the tubes on the magnet for 3 minutes.
- 9** Remove the supernatant, and replace with 500 uL of 80% EtOH. Wash for 1 min.
- 10** Remove first wash and replace with a second wash of 80% EtOH. Wash again for 1 min.
- 11** Remove second wash. Spin down the beads and any remaining EtOH, then place the tube back on the magnet and remove any remaining EtOH with a p20.
- 12** Allow the beads to dry for 2 minutes with the tubes left open.
- 13** After drying, add the desired volume of Qiagen Buffer EB or nuclease free water (generally 15 uL). Flick the tube to mix well.
- 14** Let the DNA elute off the beads for 4 minutes at room temperature.
- 15** Place the tube back on the magnet for 2 minutes, then remove the supernatant (which contains the DNA) to a new tube.
If removing the supernatant carries some of the beads over, then pipette smaller volumes separated by a minute or two that allows the magnet to pull the beads back up.