

sgRNA Cloning Protocol:

1. Anneal the two oligonucleotides together.
 - a. Resuspend the oligonucleotides to a final concentration of 100 μ M.
 - b. Mix the following:

<u>Reagent</u>	<u>Vol.</u>
100 μ M Forward Oligo	4.5 uL
100 μ M Reverse Oligo	4.5 uL
10X NEB Buffer 2	1 uL
<u>Total</u>	<u>10 uL</u>

- c. Boil 1L of water and then insert the tube containing the above mixture. Boil the mixture for 4 minutes then turn off the heat.
 - d. Allow the water to cool very slowly to room temperature (i.e. leave it on your bench at room temperature for several hours).
2. Digest the MLM3636 (Addgene #43860) plasmid with BsmBI. Gel purify the vector backbone for cloning.
3. Dilute the annealed sgRNA 1:100 and use 1 uL of the dilution for each 50 ng of vector. This produces an approximately 3:1 sgRNA:plasmid molar ratio when ligating.
4. Ligate with NEB Quick Ligase.
5. Transform into NEB 10 β competent cells.
6. Isolate the plasmid and confirm correct cloning by Sanger sequencing. We used the recommended MLM3636 sequencing primer: 5' – CAGGGTTATTGTCTCATGAGCGG – 3'. For Sanger sequencing, we submit our samples to GeneWiz.