

LeChatelier's Principle: Iron(III) Thiocyanate Equilibria


Looks like Tang, tastes like going to the hospital to get your stomach pumped

Chemicals and Equipment Needed

- LeChatelier's Principle Kit – **O2**
 - Dropper bottle of 0.10 M $\text{Fe}(\text{NO}_3)_3$
 - Dropper bottle of 0.10 M KSCN
 - Dropper bottle of 0.10 M AgNO_3
 - Small vial of NaF (s)
 - Small vial of Na_2SO_3 (s)
- ~500 mL 0.001 M $\text{Fe}(\text{NO}_3)_3$ – **H4**
 - If low, can be made from stock solution (0.1 M, also **H4**)
- ~500 mL 0.001 M KSCN – **H4**
 - If low, can be made from stock solution (0.1 M, also **H4**)
- 7-300 mL tall beakers – **Q2**
- 6 stirring rods – **U1**
- 2 microspatulas – **U1**
- Light box – **A4**

Preparation

- **MAKE SURE TO USE THE CORRECT CONCENTRATIONS!**
- Set up 7 beakers in a row.
- Add 75 mL 0.001 M $\text{Fe}(\text{NO}_3)_3$ to the 1st beaker (colorless), cover with labeled petri dish
 - Just use the marking on the beakers, no need to use graduated cylinders
- Add 75 mL 0.001 M KSCN to the 2nd beaker (colorless), cover with labeled petri dish
- Add 75 mL 0.001 M $\text{Fe}(\text{NO}_3)_3$ AND 75 mL 0.001 M KSCN to the last five beakers, producing the orange colored species FeSCN^{2+} . Cover with petri dishes labeled **FeSCN^{2+}** .
 - The solution should look like Tang or Kool-Aid. If it is dark or brown, you used the wrong concentration
- Set the dropper bottles, vials of reagents, stirring rods, and microspatulas on top of the light box according to the diagram below:




Presentation

Before each beaker, ask the audience to predict how equilibrium will shift

Beaker 1- reference color

- Add the 75 mL 0.001 M KSCN to the beaker containing only 0.001 M Fe(NO₃)₃ to show formation of the colored FeSCN²⁺ complex ion from two colorless solutions.



Beaker 2- Add a squirt or two of 0.10 M Fe(NO₃)₃


- **Observation:** The color intensifies
- **Conclusion:** There must have been free SCN⁻ available to react with the added Fe³⁺

Beaker 3- Add a squirt or two of 0.10 M KSCN

- **Observation:** The color intensifies
- **Conclusion:** There must have been free Fe³⁺ available to react with the added SCN⁻


Beaker 4- Add a small amount of solid Na₂SO₃

- **Observation:** The color fades
- **Conclusion:** The equilibrium has shifted to the left, due to removal of Fe³⁺ by reduction with SO₃²⁻




Beaker 5- Add a small amount of solid NaF

- **Observation:** The color fades
- **Conclusion:** The equilibrium has shifted to the left, due to removal of Fe³⁺ by formation of a complex with F⁻



Beaker 6- Add a squirt or two of 0.10 M AgNO₃

- **Observation:** The color fades, and the solution turns cloudy
- **Conclusion:** The equilibrium has shifted to the left, due to removal of SCN⁻ by precipitation



Clean-Up

- All the solutions can go down the sink except the one containing AgSCN (the cloudy one), which should go in the WWC

NOTES:

- To make the stock solutions:
 - 0.100 M Fe(NO₃)₃
 - 40.402g Fe(NO₃)₃•9 H₂O in 500 mL d-H₂O, add 63mL conc. HNO₃, dilute to 1L
 - 0.100 M KSCN
 - 9.718g KSCN in 500 mL d-H₂O, then dilute to 1L
- To make the 0.0010 M solutions:
 - Take 10 mL of the appropriate solution, dilute to 1L

LeChatelier's Principle: Instructor Notes

Beaker 1- reference color

- Add the 75 mL 0.001 M KSCN to the beaker containing only 0.001 M $\text{Fe}(\text{NO}_3)_3$ to show formation of the colored FeSCN^{2+} complex ion from two colorless solutions.



Beaker 2- Add a squirt or two of 0.10 M $\text{Fe}(\text{NO}_3)_3$


- **Observation:** The color intensifies
- **Conclusion:** There must have been free SCN^- available to react with the added Fe^{3+}

Beaker 3- Add a squirt or two of 0.10 M KSCN

- **Observation:** The color intensifies
- **Conclusion:** There must have been free Fe^{3+} available to react with the added SCN^-


Beaker 4- Add a small amount of solid Na_2SO_3

- **Observation:** The color fades
- **Conclusion:** The equilibrium has shifted to the left, due to removal of Fe^{3+} by reduction with SO_3^{2-}



Beaker 5- Add a small amount of solid NaF

- **Observation:** The color fades
- **Conclusion:** The equilibrium has shifted to the left, due to removal of Fe^{3+} by formation of a complex with F^-



Beaker 6- Add a squirt or two of 0.10 M AgNO_3

- **Observation:** The color fades, and the solution turns cloudy
- **Conclusion:** The equilibrium has shifted to the left, due to removal of SCN^- by precipitation

