

Barcoding for CyTOF

Written by: Laura Johnston

Date: 2/9/2018

Revised:

- **Standard protocol (Dead cell barcoding)**
 - Simplified procedure
 - Stain single cell suspension with cisplatin
 - Fix/perm
 - Stain with Pd barcoding reagents
 - Up to 20 samples
 - Combine barcoded samples into one tube
 - Fc block/stain with antibodies
 - Stain with Ir-intercalator
 - Benefits
 - Faster sample acquisition
 - Less run time between sample 1 and sample 100, minimizes machine variability between samples
 - Can collect data on samples with low cell numbers
 - Saves on reagents for staining
 - Simplifies staining step to one tube
 - Precautions
 - Cells are fixed and permeabilized before staining with antibodies, which could change the staining pattern of some markers
 - It is possible to work out which antibodies and clones work on fix/permed cells using standard flow cytometry
 - After titrating antibodies on 3×10^6 cells, the titration should be tested on the number of cells in the barcoded sample (ex: 20 samples \times 3×10^6 cells = 60×10^6 cells)
- **Modified standard protocol**
 - Simplified procedure
 - Stain single cell suspension with cisplatin
 - Fc block/stain with antibodies
 - Fix/perm
 - Stain with Pd barcoding reagents
 - Up to 20 samples
 - Combine barcoded samples into one tube
 - Stain with Ir-intercalator
 - Benefits
 - Faster sample acquisition
 - Less run time between sample 1 and sample 100, minimizes machine variability between samples
 - Can collect data on samples with low cell numbers
 - Allows for staining of cells before fix/perm if the panel contains many markers that are susceptible to fix/perm
 - Precautions
 - If you are staining 20 samples, they must be stained and barcoded individually before combining – increases the possibility of error (unevenly stained samples, mixing up tubes, etc.)
 - Does not save reagents
- **Live cell barcoding**
 - Simplified procedure
 - Stain single cell suspension with cisplatin
 - Stain with CD45 antibodies – one per sample
 - Limited by number of CD45 antibodies with distinct metal tags

- Combine CD45-tagged samples into one tube
- Fc block/stain with antibodies
- Fix/perm
- Stain with Ir-intercalator
- Benefits
 - Faster sample acquisition
 - Less run time between sample 1 and sample 100, minimizes machine variability between samples
 - Can collect data on samples with low cell numbers
 - Allows for staining of cells before fix/perm if the panel contains many markers that are susceptible to fix/perm
- Precautions
 - Limited by number of CD45 antibodies with distinct metal tags
 - Each CD45 antibody takes away a channel that could be used for your markers
 - Pd barcoding metals can be conjugated to CD45 antibodies to avoid using up channels for markers, however special custom conjugations are required. Extra troubleshooting will be required.