

Basic Flow Cytometry Staining Protocol

Written by: Laura Johnston

Date: 12/20/17

Modified: 5/15/19

Materials

- Cells to be stained
- FACS buffer – can be cell dependent
 - 1X Ca/Mg²⁺ free PBS
 - Either FBS (1-10%) or BSA (0.1-1%)
 - *Optional*: EDTA (0.5-5mM) if you have very sticky cells
 - Commonly used: 1% FBS in 1X PBS (consider filtering FBS to remove debris)
- Antibodies
- Live/dead stain (this protocol is for unfixable stains such as DAPI, PI, and the Sytox dyes)
- 5 mL flow tubes: Falcon #352008

Procedure

1. Prepare single cell suspension
2. lyse RBCs (optional, if needed)
3. Count cells
4. Adjust cell concentration to 0.1-5 million cells per 50 μ L FACS buffer
 - a. Number of cells needed depends on your experiment. General starting point is often 1 million cells.
5. Pipette 50 μ L cells into Eppendorf or 5 mL flow tubes for staining
 - a. Need one tube per sample and single stain controls
6. Add Fc block
 - a. 1-5 μ L per sample following manufacturer's protocol
 - i. Anytime you need to pipette 1 μ L it's easier to make a 1:10 dilution in FACS buffer and add 10 μ L to each sample
7. Make sure samples are mixed by gently vortexing or pipetting up and down
8. Incubate at 4°C for 10-15 minutes
9. Add antibody cocktail – 50 μ L total volume per tube
 - a. For example, if you are adding 2 μ L of antibody per tube, you would first combine 2 μ L antibody and 48 μ L FACS buffer, then add the 50 μ L of antibody cocktail to the 50 μ L of cells for a total staining volume of 100 μ L
 - b. For samples: make master mix (example below of 10X master mix) of all antibodies so that appropriate amount is in 50 μ L (when added to cells, final staining volume is 100-110 μ L, calculate for 100 μ L final volume)

| | 1X | 10X |
|--------------|----|-----|
| FITC | 1 | 10 |
| APC | 5 | 50 |
| PE | 2 | 20 |
| FACS Buffer | 42 | 420 |
| TOTAL | 50 | 500 |

- i.
 - ii. Always add an extra 1-2 to your master mix to account for pipetting error, so a 10x master mix will be good for 8-9 samples
10. Mix well and incubate for 30 minutes in the dark at 4°C (cover tubes in foil if needed)
 11. Wash cells by adding FACS buffer
 - a. 1 mL if staining in Eppendorf tubes
 - b. 2 mL if staining in flow tubes
 12. Centrifuge

13. Aspirate, pipette, or dump tubes to remove supernatant
 - a. Be careful to not disturb the pellet!
14. Resuspend in FACS buffer
 - a. Volume depends on the number of cells
 - b. Generally no less than 350 μL
15. If cells were stained in Eppendorf tubes, transfer to 5 mL FACS tubes
16. A few minutes before running samples on the cytometer, add live/dead stain (see manufacturers protocol)