Contents lists available at ScienceDirect

ELSEVIER

Journal of Pharmaceutical and Biomedical Analysis

journal homepage: www.elsevier.com/locate/jpba

Short communication

Liquid chromatography–mass spectrometry analysis of diethylcarbamazine in human plasma for clinical pharmacokinetic studies

Mark S. Schmidt^a, Christopher L. King^b, Edward K. Thomsen^c, Peter M. Siba^d, Nelly Sanuku^d, Lawrence Fleckenstein^a,*

^a College of Pharmacy, University of Iowa, Iowa City, IA 52242, USA

^b School of Medicine, Case Western Reserve University, Cleveland, OH 44106, USA

^c Liverpool School of Tropical Medicine, Liverpool L3 5QA, UK

^d Papua New Guinea Institute for Medical Research, Papua New Guinea

ARTICLE INFO

Article history: Received 14 March 2014 Received in revised form 9 May 2014 Accepted 13 May 2014 Available online 22 May 2014

Keywords: Diethylcarbamazine Liquid chromatography-mass spectrometry Plasma Antifilarial

ABSTRACT

A sensitive and selective liquid chromatographic method using mass spectrometric detection was developed for the determination of diethylcarbamazine (DEC) in human plasma. DEC and its stable isotope internal standard d₃-DEC were extracted from 0.25 mL of human plasma using solid phase extraction. Chromatography was performed using a Phenomenex Synergi 4 μ Fusion-RP column (2 mm × 250 mm) with gradient elution. The retention time was approximately 4.8 min. The assay was linear from 4 to 2200 ng/mL. Analysis of quality control samples at 12, 300, and 1700 ng/mL (*N*=15) had interday coefficients of variation of 8.4%, 5.4%, and 6.2%, respectively (*N*=15). Interday bias results were -2.2%, 6.0%, and 0.8%, respectively. Recovery of DEC from plasma ranged from 84.2% to 90.1%. The method was successfully applied to clinical samples from patients with lymphatic filariasis from a drug-drug interaction study between DEC and albendazole and/or ivermectin.

© 2014 Elsevier B.V. All rights reserved.

1. Introduction

Diethylcarbamazine (DEC, Fig. 1) is an anti-parasitic agent utilized in the treatment of lymphatic filariasis (LF). Diethylcarbamazine not only kills the microfilaria, but is the most potent available drug against adult *Wuchereria bancrofti, Brugia malayi* and *Brugia timori*, the three lymphatic dwelling parasites that infect humans [1]. In areas where onchocerciasis and loaisis are not coendemic with LF (endemic areas outside of sub-Saharan Africa), DEC combined with albendazole are the primary drugs used for annual mass drug administration (MDA) as the primary tool to elimination of LF in certain geographic areas [2]. The World Health Organization has set a goal to eradicate LF by 2020, using this MDA program with antifilarial medications with the aim to interrupt transmission and stop the spread of infection. The recommended

Abbreviations: DEC, diethylcarbamazine; LF, lymphatic filariasis; LLOQ, lower limit of quantitation; SPE, solid phase extraction.

* Corresponding author. Tel.: +1 319 335 8804; fax: +1 319 353 5646. *E-mail address*: l-fleckenstein@uiowa.edu (L. Fleckenstein).

http://dx.doi.org/10.1016/j.jpba.2014.05.016 0731-7085/© 2014 Elsevier B.V. All rights reserved. regimens must be administered once a year for at least 5 years. Although numerous clinical trials for the treatment of lymphatic filariasis have been conducted, pharmacokinetic studies of DEC in body fluids have been hampered by the lack of analytical methods with low quantitation limits and low required sample volumes. Chromatographic methods with UV detection have been published for determining DEC in medicated salts [3] or tablet formulations [4]. However, these UV methods had lower limits of quantitation $(LLOQ) \ge 100 \text{ ng/mL}$. Several gas chromatographic methods have been used to analyze DEC in plasma [5–10]. The most sensitive gas chromatographic method has a LLOQ of 10 ng/mL with a required sample volume of 1 mL of plasma using a NPD detector [9]. However, the range of the method was limited, with an upper limit of quantitation of 200 ng/mL, which necessitates many re-runs for over-range samples. A modification of the method used a FID detector and 0.5 mL of plasma has also been published; the range of the calibration curve in this approach ran from 100 to 2000 ng/mL [10].

To better support pharmacokinetic analyses for drug-drug interaction studies, a method is needed that reduces the amount of plasma needed and achieves a sensitivity of ≤ 10 ng/mL. DEC







Fig. 1. Structure of diethylcarbamazine.

is a weak base that would be expected to readily form a protonated species under typical LC–MS conditions, the detection of which should be quite sensitive relative to UV detection. Previously, gas chromatographic methods for DEC were developed because of the lack of suitable stationary phases to retain small polar compounds on liquid chromatography columns [9,10]. However, the number and quality of available stationary phases has rapidly evolved since that time, making liquid chromatography feasible for bioanalytical work with DEC. For the present work, a column/mobile phase combination was readily found with a retention factor k' of 1.5, which provided adequate resolution from endogenous substances not removed in the extraction from plasma.

2. Experimental

2.1. Solvents and chemicals

Methanol and acetonitrile were Optima grade purchased from Fisher Scientific (Fair Lawn, NJ, USA). Dibasic sodium phosphate was ACS grade purchased from Fisher Scientific. Tribasic sodium phosphate (ACS grade), formic acid (\geq 88%) and diethylcarbamazine citrate (Vetranal analytical standard) were purchased from Sigma–Aldrich (Saint Louis, MO, USA). The internal standard d₃diethylcarbamazine citrate was purchased from Toronto Research Chemicals (Toronto, Ontario, Canada). Phosphate buffer, pH 10.3 was prepared from the sodium phosphates. A solution of 0.1% formic acid in methanol was prepared to elute the solid phase extraction (SPE) cartridge. Ultrapure analytical grade Type I water was produced by a MilliQ Plus water system (Millipore Corporation, Bedford, MA, USA).

2.2. Blank plasma, plasma standards and controls, patient samples

The citrate salts of DEC were used to make stock solutions. DEC citrate and d₃-DEC citrate were weighed on a Mettler-Toledo AG104 analytical balance (Mettler-Toledo Inc., Hightstown, NJ, USA). The weights were corrected for the salt to determine the amount of free base. The weighed amount was dissolved in methanol in a volumetric flask to make 1.0 mg/mL stock solutions, which were stored at -20 °C. DEC stock solution was prepared in duplicate. Working solutions of DEC were prepared daily at 10.0, 1.0, 0.10, and 0.010 µg/mL by serial dilution of each stock solution with water. A $0.8 \,\mu g/mL$ working solution of d₃-DEC was prepared by diluting the stock solution with water in a volumetric flask, which was stored at -20 °C and used for 1 week. Calibration standards were prepared by adding between 10 and 75 µL of the appropriate working solution to a 1.5 mL microcentrifuge tube containing 20 µL of internal standard, 250 µL of plasma, and 700 µL of phosphate buffer. Control samples were prepared similarly using the duplicate set of DEC working solutions. The working curve consisted of samples containing 4, 8, 16, 64, 256, 520, 1400, and 2200 ng DEC/mL plasma. Controls were 12, 300, and 1700 ng DEC/mL plasma. Patient samples, blanks, and blank zeros were prepared by adding 20 µL of internal standard and 700 μ L of phosphate buffer to 250 μ L of patient plasma. All

solutions were vortex mixed and subjected to SPE prior to instrumental analysis. The lowest calibrant of the calibration curve was chosen because it was deemed more than sufficient for pharmacokinetic analysis and to minimize difficulties in achieving the regulatory guidelines for precision. The signal-to-noise ratio was >50 at the LLOQ, leaving some capacity to lower either the quantitation limit or the volume of plasma required for analysis if needed for specialized studies in young children or critically-ill patients.

2.3. Solid phase extraction

Agilent C18 (50 mg/1 mL) SPE Cartridges (Agilent Corp., Santa Clara, CA, USA) were used for sample preparation. A Cerex positive pressure SPE processor was used with nitrogen to modulate flow. All flow rates were approximately 1 mL/min. The cartridges were conditioned with 1 mL of methanol followed by 1 mL of water. After loading the prepared sample, the cartridge was washed with 2 mL \times 1 mL of water and eluted with 2 mL \times 1 mL of 1% formic acid in methanol. Solvent was removed under flowing nitrogen at 30 °C. The residue was reconstituted in 100 µL of 5% acetonitrile in water and transferred to an autosampler vial with glass insert for analysis. The recovery of DEC from plasma was determined by spiking triplicate plasma samples with DEC at the three control concentrations and comparing the average peak area ratio of the spiked sample to the average peak area ratio of triplicate plasma samples with analyte added post extraction, followed by drying and reconstitution. Internal standard was added post extraction in all recovery samples.

2.4. Instrumentation

The instrumentation system consisted of a Shimadzu 2010A LC–MS platform in APCI positive mode operating under LCMSSolution (Version 3.01) software (Shimadzu, Columbia, MD, USA). The analytical column was a Phenomenex Fusion-RP (4 μ m, 80 Å, 2.0 mm × 250 mm) preceded by a Phenomenex Fusion-RP Security-Guard guard (2.0 mm × 4 mm) column. Separation conditions were: sample temperature, 23 (±3)°C; column temperature, 23 (±3)°C; sample injection volume, 25 μ L. Solvent A (46%) was 0.05% formic acid; solvent B was acetonitrile with 0.05% formic acid. The analysis was run on a gradient between 5% and 57% B over 7.5 min at a flow rate of 0.25 mL/min. The total run time for a LC–MS analysis was 10 min.

The mass spectrometer was tuned using a polyethylene glycol solution following the manufacturer's protocol; the interface, CDL, and Q-array voltages optimized in the tune were used for the analysis. The scan interval was 1.0 s, microscan 0.15 amu, the APCI probe temperature was 375 °C, and the CDL (curved desolvation line) and block temperatures 200 °C. Nitrogen flow through the probe was 2.5 L/min. DEC and d₃-DEC yielded [M+1] ions with no significant adduct formation or in-probe fragmentation. The use of the stable isotope internal standard greatly lessens the likelihood of significant matrix effects.

2.5. Calculations and precision

The DEC concentrations were calculated from the peak-area ratio of DEC to d₃-DEC for standards, controls, and samples. The linear least squares equation was calculated with $1/\text{Conc}^2$ weighting. Coefficients of determination (r^2) were 0.993 or better. During validation, series of plasma controls spiked with DEC at 4, 12, 300, and 1700 ng/mL were analyzed for accuracy and precision in pentuplicate on three days (Table 1). The lower limit of quantitation was 4 ng/mL; S/N at was >50 at the LLOQ. During sample analysis,

Table 1	
Back-calculated concentrations of standard curves for DEC in human plass	ma.ª

	4	8	16	64	256	520	1400	2200
DEC concentration (ng/mL)								
Mean	3.99	8.10	15.9	62.5	257.0	538.5	1443.6	2154.9
SD	0.20	0.56	1.9	2.9	24.7	40.8	77.7	112.6
% CV	5.1	7.0	12.2	4.6	9.6	7.6	5.4	5.2
% bias	-0.3	1.2	-0.4	-2.4	0.4	3.6	3.1	-2.1

^a SD, standard deviation; CV, coefficient of variation; bias, deviation of mean from nominal value.

controls at 12, 300, and 1700 ng/mL were run in duplicate with each batch.

2.6. Clinical application

This assay was used for a pharmacokinetic drug-drug interaction study of DEC in 24 adult patients infected with lymphatic filariasis with 12 subjects in each of two arms of the study. In one arm of the study patients received a single 6 mg/kg dose of DEC citrate along with 200 mg of albendazole and in the other arm 200 μ g/kg ivermectin plus 6 mg/kg dose of DEC citrate and 200 mg of albendazole. A pre-dose sample was collected for each patient, as well as eleven time points from 1 to 72 h post-dose.

3. Results and discussion

Back calculation of the calibration curve during validation from 4 to 2200 ng/mL produced coefficients of variation between 4.6% and 12.2% and deviations from the nominal concentration (bias) of -2.4 to 3.6% (Table 1). Intraday and interday coefficients of variation for controls during validation were <12.0% at the LLOQ of 4 ng/mL and \leq 9.6% at 12, 300, or 1700 ng/mL; the intraday and interday bias was \leq 13.6% at the LLOQ and \leq 12.9% at 12, 300, or 1700 ng/mL (Table 2). Autosampler stability was assessed by analyzing a nominally 20 ng/mL sample three days consecutively, calculating the concentration from each day's calibration curve. After 30 h, the difference in concentration from the original was -11.4%. Freezer stability studies at -20°C were previously done and where plasma samples were demonstrated to be stable at

Table 2

Precision and accuracy for DEC in human plasma during validation.^a

Table 3

Stability of DEC in human plasma. Data presented as % remaining analyte in stored samples compared to initial conditions (average ± SD).

Analyte nominal concentration (ng/mL)	% remaining in plasma					
Freeze and thaw stability $N = 3^{a}$						
12	92.5 ± 4.2					
300	102.6 ± 5.8					
1700	88.4 ± 2.3					
Short-term bench stability (2 h at room temp.) N = 5						
12	89.6 ± 7.6					
300	109.3 ± 4.6					
1700	100.7 ± 7.1					
Short-term bench stability (4 h at room temp.) N = 5						
12	95.2 ± 12.1					
300	108.7 ± 3.3					
1700	100.9 ± 2.6					
Short-term bench stability (24 h at room temp.) N = 5						
12	102.5 ± 7.6					
300	98.3 ± 2.4					
1700	91.3 ± 4.2					

^a After 3 cycles to −80 °C.

-20 °C for 12 weeks [10]. Room temperature and freeze-thaw stability studies were conducted (Table 3). DEC appears to be satisfactorily stable at room temperature for at least 24 h and three freeze-thaw cycles did not cause significant reduction in drug levels. Percent recoveries of DEC from plasma were 85.1%, 84.2% and 90.1% at concentrations of 12, 300 and 1700 ng/mL, respectively (n=3).

During sample analysis of 24 subjects, all but one sample was greater than the LLOQ through 72 h. Prior studies of DEC elimination have only extended to 48 h because of quantitation limits or plasma volumes needed. The current refinements allow for better pharmacokinetic estimates as nearly all subjects will have a complete data set and the mean $AUC_{0-\infty}$ is within 2.5% of AUC_{0-t} . Fig. 2a–d are illustrative chromatograms of DEC in a blank, a patient predose, extracted standard at the LLOQ, and a patient 72 h post-dose. Fig. 3 shows a representative plasma concentration–time curve from a representative study patient. No analytical interferences were observed from other administered study drugs (albendazole and ivermectin).

		Accuracy and precision of controls					
		Day 1 (<i>n</i> =5)	Day 2 (n=5)	Day 3 (<i>n</i> =5)	Overall $(n = 15)$		
	Mean	4.30	4.54	3.75	4.20		
4 m m / mm I	SD	0.46	0.29	0.44	0.51		
4 lig/lilL	CV (%)	10.6	6.4	11.7	12.0		
	Bias (%)	7.5	13.6	-6.2	5.0		
	Mean	12.5	11.9	10.8	11.7		
12 ng/ml	SD	1.0	0.8	0.2	1.0		
12 lig/iiiL	CV (%)	8.2	6.5	2.1	8.4		
	Bias (%)	4.0	-1.0	-9.7	-2.2		
	Mean	310.0	305.0	338.7	317.9		
200 m m/m I	SD	10.5	7.2	7.4	17.3		
300 lig/lilL	CV (%)	3.4	2.4	2.2	5.4		
	Bias (%)	3.3	1.7	12.9	6.0		
	Mean	1710.0	1737.5	1691.5	1713.0		
1700 m m/m I	SD	88.3	65.4	161.8	106.4		
1700 lig/liiL	CV (%)	5.2	3.8	9.6	6.2		
	Bias (%)	0.6	2.2	-0.5	0.8		

^a SD, standard deviation; CV, coefficient of variation; bias, deviation of mean from nominal value.



Fig. 2. Chromatograms of (a) extracted blank; (b) patient sample pre-dose; (c) extracted standard at the LLOQ; (d) patient sample 72 h post-dose. Lower trace is internal standard d₃-DEC.



Fig. 3. DEC plasma concentration vs. time in a representative patient with lymphatic filariasis who received oral DEC 6 mg/kg and albendazole 200 mg given with breakfast.

4. Conclusions

We describe a LC–MS assay using solid phase extraction to quantify DEC in human plasma samples. As compared to presently published methods, the assay reduces the needed plasma sample volume from 0.5 mL to 0.25 mL, improves the quantitation limit 25-fold, and halves the time of the chromatographic run. The assay has been validated for accuracy, precision, stability, and recovery; the calibration curve is linear over the concentration range of 4–2200 ng/mL with recovery from plasma ranging from 84% to 90%. The ability of this LC–MS assay to accurately and precisely quantitate DEC in human plasma up to 72 h following a typical 6 mg/kg oral dose increases the applicability of this assay, allowing a broader range of laboratories to perform pharmacokinetic studies without the technical challenges of gas chromatography or need for much larger plasma volumes.

Acknowledgements

This research was funded in part by the Bill & Melinda Gates Foundation (GH5342, the DOLF Project). The funders had no role in study design, data collection and analysis, decision to publish, or preparation of the manuscript. The findings and conclusions contained within are those of the authors and do not necessarily reflect positions or policies of the Bill & Melinda Gates Foundation.

References

- [1] E.A. Ottesen, Rev. Infect. Dis. 7 (1985) 341-356.
- [2] E.A. Ottesen, B.O. Duke, Bull. World Health Organ. 75 (1997) 491–503.
- [3] N. Mathew, M. Kalyanasundaram, Acta Trop. 80 (2001) 97–102.
- [4] J. Mahesh Reddy, M.R. Jeyaprakash, K. Madhuri, S.N. Meyyanathan, K. Elango, Indian J. Pharm. Sci. 73 (2011) 320–323.
- [5] J.A. Bogan, Analyst 102 (1977) 56-61.
- [6] G.D. Allen, T.M. Goodchild, B.C. Weatherley, J. Chromatogr. 164 (1979) 521–526.
- [7] S. Nene, B. Anjaneyulu, T.G. Rajagopalan, J. Chromatogr. 308 (1984) 334–340.
- [8] G. Edwards, K. Awadzi, A. Breckenridge, H.M. Gilles, M.L. Orme, S.A. Ward, Clin. Pharmacol. Ther. 30 (1981) 551–557.
- [9] S. Lee, D.A. Casteel, L. Fleckenstein, J. Chromatogr. B 704 (1997) 181–185.
- [10] J.R. Miller Jr., L. Fleckenstein, J. Pharm. Biomed. Anal. 26 (2001) 665-674.